

## Project 6070: BioID hyperlopit TMT expression

Daichi Shonai, Scott Soderling

### Report for TMT Differential Protein Analysis by LC-MS/MS (7 Nov 2022)

#### Objective:

Examine the protein compositions from different centrifuged fractions (fraction4 ~ 10) of three independent study groups, each of which has 7 fractions from both WT and transgenic BioID labeled (Tg).

#### Proteomics Facility Scientific Contributors:

Duke Proteomics Core Facility: Tricia Ho (sample preparation, data collection, data analysis, report writing), Greg Waitt (data analysis), Erik J. Soderblom (study design, data collection, data analysis, report writing)

**Sample Preparation.** The Duke Proteomics and Metabolomics Core Facility (DPMCF) received 42 samples (3 of each WT\_F4, WT\_F5,... WT\_F10 and Tg\_F4, Tg\_F5,... Tg\_F10) which were kept at -80C until processing. Samples were supplemented to 3% SDS and were spiked with 1 or 2pmol casein as in internal quality control standard. They were reduced with 10 mM dithiothreitol for 30 min at 32C and alkylated with 20 mM iodoacetamide for 45 min at room temperature, then supplemented with a final concentration of 1.2% phosphoric acid and 437  $\mu$ L of S-Trap (Protifi) binding buffer (90% MeOH/100mM TEAB) (**Table 1**). Proteins were trapped on the S-Trap micro cartridge, digested using 20 ng/ $\mu$ L sequencing grade trypsin (Promega) for 1 hr at 47C, and eluted using 50 mM TEAB, followed by 0.2% FA, and lastly using 50% ACN/0.2% FA. All samples were then lyophilized to dryness.

**TMT Labeling Procedure.** Each sample was resuspended in 85  $\mu$ L 200 mM triethylammonium bicarbonate, pH 8.0 (TEAB). 65  $\mu$ L were taken of each sample and the remaining was combined to form an SPQC pooled sample. Fresh TMT16plex reagents (Thermo, 0.5 mg for each 16-plex reagent) was resuspended in 20  $\mu$ L 100% acetonitrile (ACN) and two sets were split among the three sample groups. Samples were incubated for 1 hour at RT. After 1-hour reaction, 3.25  $\mu$ L of 5% hydroxylamine was added and incubated for 15 minutes at room temperature to quench the reaction. Sample were combined then lyophilized to dryness.

**Offline Fractionation Procedure.** Each set of samples was subjected to Pierce™ High pH Reversed-Phase Peptide Fractionation Kit (ThermoFisher Scientific, see [https://www.thermofisher.com/document-connect/document-connect.html?url=https://assets.thermofisher.com/TFS-Assets%2FManuals%2F84868\\_highph\\_rp\\_peptidefract\\_UG.pdf](https://www.thermofisher.com/document-connect/document-connect.html?url=https://assets.thermofisher.com/TFS-Assets%2FManuals%2F84868_highph_rp_peptidefract_UG.pdf)). Fractions 1 and 7; Fractions 2 and 8 were combined, resulting in a total of 6 fractions per sample set. Samples were frozen and lyophilized to dryness.

**LC-MS/MS Analysis.** Quantitative LC/MS/MS was performed on 3  $\mu$ L (25%) for each fraction, using an MClass UPLC system (Waters Corp) coupled to a Thermo Fusion Lumos high resolution accurate mass tandem mass spectrometer (Thermo) equipped with a FAIMS ion-mobility device via a nanoelectrospray ionization source. Briefly, the sample was first trapped on a Symmetry C18 20 mm x 180  $\mu$ m trapping column (5  $\mu$ L/min at 99.9/0.1 v/v water/acetonitrile), after which the analytical separation was performed using a 1.8  $\mu$ m Acquity HSS T3 C18 75  $\mu$ m x 250 mm column (Waters Corp.) with a 90-min linear gradient of 5 to 30% acetonitrile with 0.1% formic acid at a flow rate of 400 nanoliters/minute (nL/min) with a column temperature of 55C. Data collection on the Fusion Lumos

mass spectrometer was performed in a data-dependent acquisition (DDA) mode with an  $r=120,000$  (@  $m/z$  200) full MS scan from  $m/z$  375 – 1575 with a target AGC target of 100% ions was performed at CVs -40v, -60v and -80v. MS/MS scans were acquired at  $r=50,000$  with AGC target of 200% and max accumulation time of 120 ms. The total cycle time was 1s per CV voltage with total cycle times of 3 sec between like full MS scans. A 45s dynamic exclusion was employed to increase depth of coverage. The total analysis cycle time for each sample injection was approximately 2 hours.

**Quantitative Data Analysis.** Following UPLC-MS/MS analyses, data were imported into Proteome Discoverer 3.0 (Thermo Scientific Inc.) where quantitative signals correlating the TMT16 reporter ions ( $m/z$  126, 127N, 127C, etc.) were extracted. In addition to quantitative signal extraction, the MS/MS data was searched against the SwissProt *M. musculus* database (downloaded in Aug 2022), a common contaminant/spiked protein database (bovine albumin, bovine casein, yeast ADH, etc.), and an equal number of reversed-sequence “decoys” for false discovery rate determination. Sequest database search algorithm was utilized to perform the database searches. Database search parameters included fixed modification on Cys (carbamidomethyl), Lys (TMT), peptide N-termini (TMT) and variable modification on Met (oxidation). Precursor mass tolerances were 2.0 ppm and product ion mass tolerances were 0.02 da with full trypsin enzyme rules required. Peptide Validator and Protein FDR Validator nodes in Proteome Discoverer were used to annotate the data at a maximum 1% protein false discovery rate based on q-value calculations. Note that peptide homology was addressed using razor rules in which a peptide matched to multiple different proteins was exclusively assigned to the protein has more identified peptides. Protein homology was addressed by grouping proteins that had the same set of peptides to account for their identification. A master protein within a group was assigned based on % coverage.

Raw peptide intensity values are presented in **Table 2**. Prior to imputation, a filter was employed which required that there were at least 2 measurements across all 42 samples for each peptide. Missing data points were imputed using a randomized intensity within the bottom 2% of the detectable signals. All peptides belonging to the same protein were then summed to the protein level and the datasets were subjected to a sample loading normalization. Following that initial normalization, the reference channel SPQC samples were used to normalize across the three individual TMT sets. Those normalized protein level data are shown in **Table 3**. Note that the total number of identified peptides within each TMT datasets is presented in the “peptides” column (number of peptides identified in set1, set 2, set 3).

## **Results**

**Depth of Coverage of LC-MS/MS analysis.** Following database searching and peptide scoring using Proteome Discoverer validation, the data was annotated at a 1% protein false discovery rate. The overall coverage was **8,838 peptides corresponding to 747 unique proteins**.

**Measures of analytical versus biological variability.** In order to assess technical reproducibility, we calculated % coefficient of variation (%CV) for each protein or phosphopeptide across the 6 SPQC pool channels that were interspersed throughout the study (**Table 3**). The mean %CV of the SPQC pools was **43.8%**. To assess biological + processing + technical variability, %CVs were measured for each protein which averaged **27.7%**.

**Initial Statistical Analysis.** As an initial statistical analysis, we calculated fold-changes between various groups based on the protein expression values and calculated two-tailed heteroscedastic t-test on log2-transformed data for each of these comparisons (**Tables 4-11**). The column “x\_FC” in those tables indicated fold change relative to the second group in the header. For example, if there was a protein that said 1.9 for tg vs wt, then it was 1.9-fold higher in tg relative to wt. If a protein met the

criteria of being at least 1.3 fold different between the groups with a p-value < 0.05 then we indicated that in the “Stats” column. If it says “up” this indicates the peptide was higher in tg relative to wt by at least 1.3 fold with a p-value of < 0.05.

**Conclusions.** Although the number of peptides was reasonable for this study (>8,800) we unfortunately did not achieve very high protein coverage (~750 proteins). Some of this could have been lower than ideal loading on the system from the pre-enriched fractions. Note that we fractionated using our most sensitive approach (offline C18 high pH) and loaded as much (~ 40%) on the system as we could. There were a handful of proteins that appear to change localization as a function of the tg mutation for each of the unique fractions. A follow-up could be to compare those proteins across the different fractions to see if any go down in one fraction while increasing in another (indicating translocation). Also note that we chose an arbitrary fold-change cutoff of 1.3. There were other proteins passing the 0.05 pvalue threshold that were not indicated as “up” or “down” that you could consider.

**Supplemental Data.** This document (6070\_Report\_110722.pdf) and an Excel spreadsheet containing Tables 1-11 (6070\_SupplementalData\_110722.xlsx) has been uploaded to the Project 6070 folder of the Express Data Repository (<https://discovery.genome.duke.edu/express/resources/6070/>).

**Acknowledgment.** As the Fusion Lumos Tribrid mass spectrometer used to generate this data was funded in part by a S10 award, please include the following as an acknowledgement in any publication or presentation made based on this dataset:

"Research reported in this publication was supported in part by the Office of The Director, National Institutes of Health of the National Institutes of Health under Award Number S10OD024999. The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health." Also, please inform the DPMSR of any publications from this instrument as these have to be included in our annual report` to the NIH.